EN Instructions for Use
trophon® NanoNebulant®

35% hydrogen peroxide

Single use disinfectant delivery cartridge, containing multiple doses of disinfectant.
The following must be read in conjunction with your trophon® EPR, trophon®2, or trophon®3 User Manual. This cartridge is specifically designed for use with the trophon EPR, trophon2, and trophon3 to deliver the correct dosage of disinfectant to ensure high level disinfection (HLD).
Each cartridge contains 80mL of 35% w/w hydrogen peroxide solution. Use of disinfectants with differing composition or concentration to that specified may result in damage to the trophon EPR, trophon2, trophon3, probe and/or reduced effectiveness of the HLD process.

The disinfectant cartridge is intended to be used only in combination with the trophon EPR, trophon2 and trophon3 device. Once inserted into trophon EPR, trophon2 or trophon3 device, a high level instrument-grade disinfectant system is formed.

Refer to your trophon EPR, trophon2 or trophon3 User Manual to enable the AcuTrace® Technology, if available.

Intended Use
The trophon EPR, trophon2 and trophon3 are designed to provide High-Level Disinfection (HLD) of validated ultrasound probes. High-Level Disinfection is achieved by surface exposure to a controlled dose of hydrogen peroxide mist delivered to a disinfection chamber containing the ultrasound probe.

The trophon EPR, trophon2 and trophon3 system consist of a multiple-use instrument combined with a single use disinfectant "trophon NanoNebulant", delivered from a multi-dose cartridge.

1.Warnings (Refer to the Safety Data Sheet for more information)
IMPORTANT: Cartridges are pierced at the top and the bottom side upon insertion of a new cartridge into the trophon EPR, trophon2 and trophon3. A small amount of disinfectant may remain in the cartridge, even when it has been fully used. Follow the instructions carefully to avoid injury.
• Always wear disposable gloves when handling disinfectant cartridges.
• Avoid contact with eyes and skin.
• Do not manually pierce the cartridge.
• Do not attempt to open or use a damaged or distorted cartridge.
• Do not attempt to refill or reprocess the empty cartridge after it has been removed from the trophon EPR, trophon2 and trophon3. This could lead to personnel injury and/or equipment/trophon EPR, trophon2 and trophon3 malfunction or damage.
• Strong Oxidizer: Contact with other materials may cause a fire.
• Corrosive: Harmful if swallowed or inhaled. May cause eye and skin burns. Prolonged or repeated skin contact may cause dermatitis.

Any serious incidents in relation to trophon2 or trophon3 or trophon NanoNebulant should be reported to Nanosonics and / or the authority of your member state.

2.Emergency Guidelines
• **Spillage:** Always wear personal protective equipment. Turn OFF the trophon EPR, trophon2 and trophon3 at the power switch and at the mains power switch. Dilute with plenty of water and discard according to Occupational Health and Safety Guidelines for your institution.
• **Eye contact:** Rinse immediately with copious amounts of cold water, including under the eyelids for at least 15 minutes. Keep eyes open and do NOT rub eyes.
• **Skin contact:** Flush skin with copious amounts of soap and water for at least 15 minutes. Remove contaminated clothing and shoes. Consult a physician if irritation persists.
• **Ingestion:** Do NOT induce vomiting. If victim is conscious and alert, rinse mouth thoroughly with water. Never give anything by mouth to an unconscious person. Call a poison control centre.
• **Inhalation:** Remove person to fresh air. If symptoms persist, call a physician.

3.Warning Symbol Definitions

	Oxidiser		Warning		Corrosive
	Environmental Hazard		Wear Gloves		Dangerous goods packed in limited quantities
	Caution	MD	Medical Device		Single Use Only

4.Other Symbol Definitions

LOT	Batch No.	REF	Product No.		Consult Instructions for Use
AcuTrace RFID Zone	AcuTrace® RFID Zone		Expiry Date	EU REP	European Authorised Representative
	EU Importer	CE 0197	Conforming to Medical Devices Regulation (EU)2017/745 (MDR)		
	No Air Transport	UK CA 2571	Conforming to UK MDR 2002 (SI 2002 No 618)		

5.Inserting the Disinfectant Cartridge
IMPORTANT: Do NOT use a distorted cartridge. Always check the expiration date on the cartridge before use.

Cartridge door will automatically open when the cartridge needs replacing. Do NOT force the cartridge door open as this may damage and render the trophon EPR, trophon2 or trophon3 non functional.

i. If AcuTrace is enabled, scan the AcuTrace RFID Zone on the cartridge label (figure 1).
ii. Remove the cap from the cartridge
iii.Place the cartridge NECK FIRST into the holder. Do NOT force the cartridge into the holder (figure 2).
iv.When situated in place correctly, the bottom of the cartridge should be in line with the top of the holder. The door should close smoothly. Do NOT force the door to close.
v.The trophon EPR, trophon2 or trophon3 is now ready for use. Once the door is closed, it will automatically lock and cannot be reopened until the cartridge is consumed or purged.
6.Cartridge Removal
Cartridge removal is required when the cartridge has expired or is empty. This will be indicated by a screen message.
The trophon EPR, trophon2 or trophon3 will need to be purged before being transported or when a cartridge has expired to ensure all disinfectant is removed. Refer to the trophon EPR, trophon2 or trophon3 User Manual for further information on purging.

i. Follow the trophon EPR, trophon2 or trophon3 onscreen instructions. Wear gloves, the door will open automatically. When directed, the empty cartridge can be lifted out. Avoid touching the pierced areas of the cartridge.
ii. Do NOT shake or change the orientation of the cartridge. For empty cartridge disposal see section 7 below.
iii. Replace with a new cartridge. Refer to section 5 above.

7.Disposal of Used Cartridge
Follow procedures of your institution for the disposal of CORROSIVE or OXIDIZING materials or dispose of as per local environmental and government regulatory requirements.

8.Storage

• Store cartridges in ALL original packaging until use.
• Store cartridge upright as indicated on cartridge packaging.
• Cartridges may distort or explode if not properly stored.

DANGER

Harmful if inhaled. Harmful if swallowed. Causes skin irritation. Causes serious eye damage. May cause respiratory irritation. May intensify fire; oxidizer.

Avoid breathing dust/ mist/vapours/spray. Wear protective gloves/ protective clothing. IF ON SKIN: Wash with plenty of soap and water. Dispose of container in accordance with local/regional/national regulation.

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